




Anti-PSME1 antibody ab3333

4 References 1 图像

概述	
产品名称	Anti-PSME1抗体
描述	兔多克隆抗体to PSME1
宿主	Rabbit
特异性	Detects PSME1 (proteasome activator complex subunit 1).
经测试应用	适用于: WB
种属反应性	与反应: Mouse  预测可用于: Cow, Pig, Cynomolgus monkey 
免疫原	Synthetic peptide corresponding to Human PSME1 aa 5-22. Sequence:  RVQPEAQAKVDVFREDLC  (Peptide available as <a href="#">ab4949</a> ) <div> <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a></div>
阳性对照	Rat brain extract.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituents: 0.1% BSA, 99% PBS
纯度	Immunogen affinity purified
克隆	多克隆

同种型IgG

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab3333于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 28 kDa (predicted molecular weight: 29 kDa). By Western blot, this antibody detects a 28 kDa protein representing PSME1 in rat brain extract.

靶标

功能	Implicated in immunoproteasome assembly and required for efficient antigen processing. The PA28 activator complex enhances the generation of class I binding peptides by altering the cleavage pattern of the proteasome.
序列相似性	Belongs to the PA28 family.

图片



Western blot detection of PSME1 in Hepa 1 cell lysate using ab3333.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors