

Anti-PVL LukS antibody ab190473


1 References 1 图像

概述	
产品名称	Anti-PVL LukS 抗体
描述	兔多克隆抗体 to PVL LukS
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Staphylococcus aureus
免疫原	Recombinant full length protein (His-tag) corresponding to PVL LukS. (Staphylococcus aureus PVL LukS purified from E.coli).
阳性对照	WB: Purified PVL Luk-S.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium azide Constituent: 99% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG
应用	

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 0.5 µg/ml. Predicted molecular weight: 35 kDa. If probing against purified protein, at least 10 ng of pure toxin can be detected with 0.5 ug/ml of the primary antibody combined with a suitable dilution of a secondary anti-rabbit antibody

## 图片



Anti-PVL LukS antibody (ab190473) at 0.5 µg/ml + Purified PVL  
Luk-S

### Secondary

Goat anti rabbit-AP at 1/3000 dilution

**Predicted band size:** 35 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors