

# Anti-Red Blood Cells antibody ab34858

## 6 References

### 概述

产品名称	Anti-Red Blood Cells抗体
描述	兔多克隆抗体to Red Blood Cells
宿主	Rabbit
经测试应用	适用于: Agglutination
种属反应性	与反应: Human
免疫原	Tissue/ cell preparation: Human washed pooled Red Blood Cells.
常规说明	<p>ab34858 is whole antiserum. For your information, we provide you the total protein concentration which is 90 mg/ml. Since the antibody is unpurified, we can only estimate its concentration at 1–10 mg/ml.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 0.88% Sodium chloride, 0.42% Potassium phosphate</p>
纯化说明	Ab34858 was prepared from polyspecific antiserum by a delipidation and defibrination.
克隆	多克隆
同种型	IgG

### 应用

## The Abpromise guarantee **Abpromise™**承诺保证使用ab34858于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Agglutination		Use at an assay dependent concentration. An agglutination titer against serial dilutions of washed pooled cells on a titer plate of 1:64 was observed for this product.

## 靶标

**相关性** Red blood cells are the most common type of blood cell and are the vertebrate body's principal means of delivering oxygen from the lungs or gills to body tissues via the blood.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors