

# Anti-Reovirus antibody ab33986

## 1 References

### 概述

产品名称	Anti-Reovirus抗体
描述	鸡多克隆抗体to Reovirus
宿主	Chicken
特异性	Monospecific. Neutralizing indices reported with each lot. Avian Reovirus AGP ELISA Positive Avian Encephalomyelitis ELISA Negative Avian Influenza (Type A) ELISA Negative Avian Rhinotracheitis ELISA Negative Infectious Bronchitis Virus ELISA Negative Newcastle Disease ELISA Negative Infectious Bursal Disease ELISA Negative Infectious Laryngotracheitis ELISA Negative Avian Reovirus AGP Positive Avian Adenovirus AGP Negative Avian Rotavirus AGP Negative Infectious Bursal Disease AGP Negative Avian Encephalomyelitis AGP Negative Fowl Pox AGP Negative Avian Adenovirus Group II (HEV)AGP Negative Avian Influenza (Type A) AGP Negative Infectious Laryngotracheitis AGP Negative Marek's Disease AGP Negative Newcastle Disease Virus AGP Negative Reticuloendotheliosis AGP Negative Avian Reovirus MNT 10e2.4 Avian Adenovirus Group I MNT Negative Avian Leukosis (RSV A&B) MNT Negative
经测试应用	适用于: Neutralising, ELISA
种属反应性	与反应: Reovirus sp.
免疫原	Reovirus
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
纯度	Whole antiserum
克隆	多克隆

## 应用

**The Abpromise guarantee** **Abpromise™**承诺保证使用ab33986于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Neutralising		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

## 靶标

## 相关性

Reovirus is an acronym for Respiratory Enteric Orphan Virus. The reovirus is a double stranded RNA virus and is non-pathogenic naturally occurring virus, and not usually associated with any particular illness. Almost 100 per cent of the population show signs of a previous reovirus infection, generally confined to the gastrointestinal or respiratory systems in the body. Scientific studies have demonstrated that the virus specifically targets cells with a highly activated Ras signaling pathway and hence is a potent cancer killer.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors