


### Anti-RFC5 antibody ab3619

[2 References](#) [1 图像](#)

#### 概述

产品名称	Anti-RFC5抗体
描述	山羊多克隆抗体to RFC5
宿主	Goat
特异性	We have data to indicate that this antibody may not cross react with Mouse. However, this has not been conclusively tested and expression levels may vary in certain cell lines/tissues.
经测试应用	<b>适用于:</b> WB
种属反应性	<b>与反应:</b> Human <b>预测可用于:</b> Chimpanzee, Gorilla, Orangutan 
免疫原	Synthetic peptide corresponding to Human RFC5.
阳性对照	HeLa cell nuclear extract.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7 Preservative: 0.1% Sodium azide Constituents: 0.021% PBS, 1.764% Sodium citrate, 1.815% Tris
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG


**The Abpromise guarantee**      **Abpromise™承诺保证使用ab3619于以下的经测试应用**

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/5000. Detects a band of approximately 37 kDa (predicted molecular weight: 36 kDa).

**靶标**

<b>功能</b>	The elongation of primed DNA templates by DNA polymerase delta and epsilon requires the action of the accessory proteins proliferating cell nuclear antigen (PCNA) and activator 1.
<b>序列相似性</b>	Belongs to the activator 1 small subunits family.
<b>细胞定位</b>	Nucleus.

**图片**

Sample: 50µg HeLa cell nuclear extract

Anti-RFC5 (ab3619) used at 1:1000 dilution

Sample: 50µg HeLa cell nuclear extract

Anti-RFC5 (ab3619) used at 1:1000 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors