

Anti-RFX6 antibody ab103497

2 图像

概述

产品名称	Anti-RFX6抗体
描述	兔多克隆抗体to RFX6
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	<p>Recombinant full length protein within Human RFX6 aa 1-950. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please <u>contact</u> our Scientific Support team to discuss your requirements.</p> <p>Database link: <u>NP_775831.1</u></p>
阳性对照	Human kidney lysate and RFX6 transfected 293T cell line lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.40</p> <p>Constituent: 100% PBS</p>
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab103497于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 102 kDa.

靶标

相关性 RFX6 belongs to the RFX family, which encodes transcription factors that contain a highly-conserved winged helix DNA binding domain.


细胞定位 Nucleus

图片



Anti-RFX6 antibody (ab103497) at 1/500 dilution + Human kidney lysate at 50 µg

Predicted band size: 102 kDa



All lanes : Anti-RFX6 antibody (ab103497) at 1/500 dilution

Lane 1 : RFX6 transfected 293T cell line lysate

Lane 2 : Non transfected 293T cell line lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 102 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors