

Anti-RPP30 antibody [2931D5a] ab53744

2 图像

概述

产品名称	Anti-RPP30抗体[2931D5a]
描述	小鼠单克隆抗体[2931D5a] to RPP30
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Recombinant fragment (Human)
阳性对照	HeLa whole cell lysate NIH3T3 whole cell lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p>
纯度	Protein G purified
纯化说明	Purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG depleted (approximately 95%) fetal bovine serum.
克隆	单克隆
克隆编号	2931D5a
同种型	IgG1

应用


The Abpromise guarantee **Abpromise™**承诺保证使用ab53744于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/100. Detects a band of approximately 30 kDa (predicted molecular weight: 30 kDa).

靶标

图片




Anti-RPP30 antibody [2931D5a] (ab53744) at 1/100 dilution +
HeLa whole cell lysate at 25 µg

Secondary
anti Mouse IgG at 1/2500 dilution

Developed using the ECL technique.

Predicted band size: 30 kDa
Observed band size: 30 kDa



Anti-RPP30 antibody [2931D5a] (ab53744) at 1/100 dilution +
NIH3T3 whole cell lysate at 25 µg

Secondary
anti Mouse IgG at 1/2500 dilution

Developed using the ECL technique.

Predicted band size: 30 kDa
Observed band size: 30 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors