

Anti-RPS9 antibody ab74711

4 References 2 图像

概述	
产品名称	Anti-RPS9抗体
描述	兔多克隆抗体to RPS9
宿主	Rabbit
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human
免疫原	Synthetic peptide derived from an internal region of human RPS9
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, PBS
纯度	Without Mg2+ and Ca2+
克隆	Immunogen affinity purified
同种型	多克隆 IgG
应用	

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 26 kDa.
ICC/IF		1/500 - 1/1000.

靶标


序列相似性

Belongs to the ribosomal protein S4P family.
Contains 1 S4 RNA-binding domain.

细胞定位

Cytoplasm. Localized in cytoplasmic mRNP granules containing untranslated mRNAs.

图片



All lanes : Anti-RPS9 antibody (ab74711) at 1/500 dilution

Lane 1 : extracts from HepG2 cells at 20 µg

Lane 2 : extracts from COLO
cells at 20 µg

Lane 3 : extracts from HT-29 cells at 20 µg


Lane 4 : extracts from HeLa cells at 30 µg

Lane 5 : extracts from HUVEC cells at 20 µg

Lane 6 : extracts from HepG2 cells at 20 µg with immunizing
peptide at 10 µg

Predicted band size: 26 kDa

Observed band size: 26 kDa



Immunofluorescence analysis of A549 cells, using ab74711 at 1/500-1/1000 dilution in the presence (+) or absence (-) of immunising peptide.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors