

Anti-S100A9 + Calprotectin (S100A8/A9 complex) antibody [MAC387] - Low endotoxin, Azide free ab80084

4 References **1 图像**

概述

产品名称	Anti-S100A9 + Calprotectin (S100A8/A9 complex)抗体[MAC387] - Low endotoxin, Azide free
描述	小鼠单克隆抗体[MAC387] to S100A9 + Calprotectin (S100A8/A9 complex) - Low endotoxin, Azide free
宿主	Mouse
经测试应用	适用于: Flow Cyt (Intra)
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human S100A9 + Calprotectin (S100A8/A9 complex). Solubilized human monocyte
阳性对照	Flow Cyt (Intra): Human peripheral blood granulocytes.
常规说明	Endotoxin Level <0.01EU/ug. The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
无载体	是
纯度	Protein G purified
克隆	单克隆
克隆编号	MAC387
骨髓瘤	NS1
同种型	IgG1

应用


The Abpromise guarantee

Abpromise™承诺保证使用ab80084于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt (Intra)		1/50 - 1/200. Membrane permeabilisation is required for this application. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

图片



ab80085 at 1/50-1/200 dilution staining S100A9 + Calprotectin in human peripheral blood granulocytes by intracellular flow cytometry.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors