

Anti-SCGB3A1 antibody ab92350

★★★★★ [1 Abreviews](#)

概述

产品名称

Anti-SCGB3A1抗体

描述

山羊多克隆抗体to SCGB3A1


① 该产品是 **快线抗体**，通过亲和纯化得到。对抗原肽有高灵敏度及特异性(ELISA测试)。

[Read the terms of use »](#)

宿主

Goat

种属反应性

预测可用于: Human 

免疫原

Synthetic peptide:

NHLIEGSQKCVAE

, corresponding to internal sequence amino acids 69 - 81 of Human SCGB3A1 (NP_443095.2)

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式

Liquid

存放说明

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

存储溶液

pH: 7.30

Preservative: 0.02% Sodium azide

Constituents: Tris buffered saline, 0.5% BSA

纯度

Immunogen affinity purified

纯化说明

ab92350 is purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

克隆 多克隆
同种型 IgG

应用

Fast track antibodies constitute a diverse group of products that have been released to accelerate your research, but are not yet fully characterized. They have all been affinity purified and show high titre values against the immunizing peptide (by ELISA).

Fast track terms of use

应用	Ab评论	说明
WB		<p>Use at an assay dependent concentration. Predicted molecular weight: 10 kDa.</p> <p>A 1 hour primary incubation is recommended for this product. Preliminary experiments gave an approx 33kDa band in Human Breast cancer and normal Lung lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation</p>

靶标

功能 Potential growth inhibitory cytokine.
组织特异性 Highly expressed in breast tissues. Absent in breast cancer cell lines.
序列相似性 Belongs to the secretoglobin family. UGRP subfamily.
细胞定位 Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors