


Anti-SLC32A1/VGAT antibody ab42939

★★★★★ [1 Abreviews](#) [3 References](#) [1 图像](#)

概述

产品名称	Anti-SLC32A1/VGAT抗体
描述	兔多克隆抗体to SLC32A1/VGAT
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Rat 预测可用于: Mouse, Cow, Dog, Non human primates 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	Rat hippocampal lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.50 Constituents: 0.238% HEPES, 50% Glycerol, 0.87% Sodium chloride, 0.01% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee [Abpromise™](#)承诺保证使用ab42939于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	1/1000. Detects a band of approximately 53 kDa (predicted molecular weight: 58 kDa).

## 靶标

功能	Involved in the uptake of GABA and glycine into the synaptic vesicles.
组织特异性	Retina. Expressed throughout the horizontal cells or more specifically at the terminals.
序列相似性	Belongs to the amino acid/polyamine transporter 2 family.
细胞定位	Cytoplasmic vesicle membrane.

## 图片



Anti-SLC32A1/VGAT antibody (ab42939) at 1/1000 dilution + Rat hippocampal lysate at 20 µg

### Secondary

Goat-anti-rabbit HRP. at 1/15000 dilution

**Predicted band size:** 58 kDa

**Observed band size:** 53 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors