

# Anti-SLITRK2 antibody ab67304

## 1 图像

### 概述

---

<b>产品名称</b>	Anti-SLITRK2抗体
<b>描述</b>	兔多克隆抗体to SLITRK2
<b>宿主</b>	Rabbit
<b>特异性</b>	This antibody is predicted to have no cross-reactivity to other SLITRK proteins.
<b>经测试应用</b>	<b>适用于:</b> WB
<b>种属反应性</b>	<b>与反应:</b> Rat
<b>免疫原</b>	Synthetic peptide (Human) of 17 amino acids from near the C-terminus
<b>阳性对照</b>	Rat Brain Tissue Lysate
<b>常规说明</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

### 应用

---

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab67304于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 95 kDa (predicted molecular weight: 95 kDa).

## 靶标

功能	Suppresses neurite outgrowth.
组织特异性	Expressed predominantly in the cerebral cortex of the brain but also at low levels in the spinal cord and medulla. Also expressed in some astrocytic brain tumors such as astrocytomas, oligodendrogliomas, glioblastomas, gangliogliomas and primitive neuroectodermal tumors.
序列相似性	Belongs to the SLITRK family. Contains 12 LRR (leucine-rich) repeats. Contains 2 LRRCT domains. Contains 1 LRRNT domain.
细胞定位	Membrane.

## 图片



**Lane 1 :** Anti-SLITRK2 antibody (ab67304) at 1 µg/ml

**Lane 2 :** Anti-SLITRK2 antibody (ab67304) at 2 µg/ml

**All lanes :** Rat brain tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 95 kDa

**Observed band size:** 95 kDa

**Additional bands at:** 70 kDa, 90 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors