

Anti-SLP76 antibody ab17029

★★★★★ [1 Abreviews](#) [1 References](#) [1 图像](#)

概述

产品名称	Anti-SLP76抗体
描述	山羊多克隆抗体to SLP76
宿主	Goat
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide: ALRNVPRSEV , corresponding to amino acids 2-12 of Human SLP76. Run BLAST with Run BLAST with
阳性对照	Recombinant Human SLP76 protein (ab114697) can be used as a positive control in WB.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: Tris buffered saline, 0.5% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab17029于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	Use a concentration of 0.5 - 2 µg/ml. Detects a band of approximately 75 kDa (predicted molecular weight: 61 kDa). 1 hour primary incubation is recommended for this product.

靶标

功能	Involved in T-cell antigen receptor mediated signaling.
组织特异性	Highly expressed in spleen, thymus and peripheral blood leukocytes. Highly expressed also in T-cell and monocytic cell lines, expressed at lower level in B-cell lines. Not detected in fibroblast or neuroblastoma cell lines.
序列相似性	Contains 1 SAM (sterile alpha motif) domain. Contains 1 SH2 domain.
结构域	The SH2 domain mediates interaction with SHB.
翻译后修饰	Phosphorylated after T-cell receptor activation by ZAP-70.
细胞定位	Cytoplasm.

图片



Western blot - Anti-SLP76 antibody (ab17029)

Anti-SLP76 antibody (ab17029) at 0.3 µg/ml + Jurkat lysate (protein in RIPA buffer) at 35 µg

Predicted band size: 61 kDa
Observed band size: 80 kDa

Primary incubation was 1 hour. Detected by chemiluminescence.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors