

Anti-SLPI antibody [31] ab17157

★★★★★ [1 Abreviews](#) [6 References](#) [1 图像](#)

概述

产品名称	Anti-SLPI抗体[31]
描述	小鼠单克隆抗体[31] to SLPI
宿主	Mouse
特异性	The reactivity of the antibody with the isolated domains was evaluated using isolated domains of SLPI obtained by cleavage using partial acidic hydrolysis. Therefore, it is possible that the antibody will recognize other SLPI cleavage products generated by e.g. (microbial) proteases.
经测试应用	适用于: IHC-Fr
种属反应性	与反应: Human
免疫原	Full length native protein (purified) corresponding to Human SLPI. SLPI purified from sputum
阳性对照	PMN lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
克隆	单克隆
克隆编号	31
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab17157于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		1/50.

靶标

功能	Acid-stable proteinase inhibitor with strong affinities for trypsin, chymotrypsin, elastase, and cathepsin G. May prevent elastase-mediated damage to oral and possibly other mucosal tissues.
组织特异性	Mucous fluids.
序列相似性	Contains 2 WAP domains.
细胞定位	Secreted.

图片



Immunohistochemical analysis on human tonsil frozen sections with ab17157 at a 1:50 dilution.

Immunohistochemistry (Frozen sections) - Anti-SLPI antibody [31] (ab17157)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors