


Anti-SOAT 1/ACAT1 antibody ab39327

15 References 2 图像

概述	
产品名称	Anti-SOAT 1/ACAT1抗体
描述	兔多克隆抗体to SOAT 1/ACAT1
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Rat, Human 预测可用于: Cynomolgus monkey 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 50% Glycerol, 0.1% BSA
纯度	Affinity purified
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee Abpromise™承诺保证使用ab39327于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 47 kDa.
IHC-P		Use at an assay dependent concentration.

靶标

功能	Catalyzes the formation of fatty acid-cholesterol esters. Plays a role in lipoprotein assembly and dietary cholesterol absorption. In addition to its acyltransferase activity, it may act as a ligase.
序列相似性	Belongs to the membrane-bound acyltransferase family. Sterol o-acyltransferase subfamily.
细胞定位	Endoplasmic reticulum membrane.

图片




All lanes : Anti-SOAT 1/ACAT1 antibody (ab39327)

Lane 1 : Rat Liver (30µg)
Lane 2 : HepG2 cell lysate (30µg)

Predicted band size: 47 kDa
Observed band size: 50,65 kDa

The 65 kDa band detected is the native form of ACAT1.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) of human adrenal tissue using ab39327 at 1/200 after heat induced antigen retrieval in pH 6.0 citrate buffer. DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SOAT 1/ACAT1 antibody (ab39327)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors