

Anti-SOX4 antibody [154C4a] ab70598

★★★★★ [2 Abreviews](#) [12 References](#) [1 图像](#)

概述

产品名称	Anti-SOX4抗体[154C4a]
描述	小鼠单克隆抗体[154C4a] to SOX4
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Recombinant fragment
免疫原	Recombinant fragment of human SOX4 from an internal sequence.
阳性对照	Recombinant human SOX4.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.1136% Dibasic monohydrogen sodium phosphate, 1% BSA, 0.03% Potassium phosphate, 0.812% Sodium chloride, 0.1312% Sodium phosphate, 0.0225% Potassium chloride, 0.0204% Monobasic dihydrogen potassium phosphate
纯度	Protein G purified
纯化说明	ab70598 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum and filtered through a 0.22µm membrane.
克隆	单克隆
克隆编号	154C4a

同种型

IgG2a

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab70598于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★☆☆☆☆ (2)	Use at an assay dependent concentration. Detects a band of approximately 38 kDa (predicted molecular weight: 47 kDa).

靶标

功能	Transcriptional activator that binds with high affinity to the T-cell enhancer motif 5'-AACAAAG-3' motif.
组织特异性	Testis, brain, and heart.
序列相似性	Contains 1 HMG box DNA-binding domain.
细胞定位	Nucleus.

图片



Anti-SOX4 antibody [154C4a] (ab70598) + immunising recombinant protein

Predicted band size: 47 kDa

Observed band size: 38 kDa

The molecular weight of the band on the western blot does not correspond to the molecular weight of the natural protein because only a fragment of the protein was used.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors