


### Anti-SOX5 antibody ab94396

★★★★★ [5 Abreviews](#) [29 References](#) [1 图像](#)

#### 概述

产品名称	Anti-SOX5抗体
描述	兔多克隆抗体to SOX5
宿主	Rabbit
特异性	Replenishment batches of ab94396 are tested in WB. Previous batches were additionally validated in ChIP and IHC-Fr. These applications are still expected to work and are covered by our Abpromise guarantee.
经测试应用	<b>适用于:</b> IHC-Fr, ChIP, WB
种属反应性	<b>与反应:</b> Mouse, Human <b>预测可用于:</b> Chicken, Cow, Pig, Zebrafish 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituent: PBS  Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.
纯度	Immunogen affinity purified

克隆

多克隆

同种型

IgG

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab94396于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr	★★★★★ (1)	Use at an assay dependent concentration.
ChIP		Use at an assay dependent concentration. PubMed: 22547066
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 84 kDa (predicted molecular weight: 84 kDa).

靶标

功能

BindsspecificallytotheDNAsequence5'-AACAAAT-3'.ActivatestranscriptionofCOL2A1andAGC1invitro.


序列相似性

Contains1HMGboxDNA-bindingdomain.

细胞定位

Nucleus.

图片



Anti-SOX5 antibody (ab94396) at 1 µg/ml + Human brain tissue lysate - total protein ([ab29466](#)) at 10 µg

**Secondary**  
Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ([ab97080](#)) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 84 kDa  
**Observed band size:** 84 kDa  
**Additional bands at:** 24 kDa, 28 kDa, 55 kDa. We are unsure as to the identity of these extra bands.

**Exposure time:** 5 minutes

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors