

Anti-SPINK1/P12 antibody [SP166] ab183347

重组 RabMAb

3 图像

概述

| | |
|-------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 产品名称 | Anti-SPINK1/P12抗体[SP166] |
| 描述 | 兔单克隆抗体[SP166] to SPINK1/P12 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, IHC-P |
| 种属反应性 | 与反应: Human |
| 免疫原 | Recombinant full length protein corresponding to Human SPINK1/P12 aa 1 to the C-terminus. Without signal peptide. Database link: P00995 |
| 阳性对照 | Human prostate adenocarcinoma tissue; Human pancreas lysate |
| 常规说明 | Cow: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information. This product is a recombinant monoclonal antibody, which offers several advantages including: - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . This product is FOR RESEARCH USE ONLY. For commercial use, please contact partnerships@abcam.com. |

性能

| | |
|------|-----------------------------------------------------------------------------------------------------------------------------------|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 7.60 Preservative: 0.1% Sodium azide Constituents: PBS, 1% BSA |
| 纯度 | Protein A/G purified |
| 纯化说明 | Purified from TCS by protein A/G. |

| | |
|------|-------|
| 克隆 | 单克隆 |
| 克隆编号 | SP166 |
| 同种型 | IgG |

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab183347于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|-------|------|---------------------------------------------------------------------------------------------------------------------------|
| WB | | 1/100. Predicted molecular weight: 8 kDa. |
| IHC-P | | 1/100. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. |

靶标


| | |
|-------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 功能 | Serine protease inhibitor which exhibits anti-trypsin activity (PubMed:7142173). In the pancreas, protects against trypsin-catalyzed premature activation of zymogens. In the male reproductive tract, binds to sperm heads where it modulates sperm capacitance by inhibiting calcium uptake and nitrogen oxide (NO) production. |
| 疾病相关 | Pancreatitis, hereditary Tropical calcific pancreatitis |
| 序列相似性 | Contains 1 Kazal-like domain. |
| 细胞定位 | Secreted. |

图片



Anti-SPINK1/P12 antibody [SP166] (ab183347) at 1/100 dilution + Human pancreas lysate

Predicted band size: 8 kDa



Immunohistochemical analysis of paraffin embedded Human prostate adenocarcinoma tissue labeling SPINK1/P12 with ab183347 at 1/100 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SPINK1/P12 antibody [SP166] (ab183347)

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-SPINK1/P12 antibody [SP166] (ab183347)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors