

Anti-STK17A antibody ab97530

[2 References](#) [2 图像](#)

概述

产品名称	Anti-STK17A抗体
描述	兔多克隆抗体to STK17A
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Cow 
免疫原	Recombinant fragment corresponding to Human STK17A aa 1-201. (AAH47696)
阳性对照	293T, A431, H1299, HeLaS3, Molt-4 and Raji cell lysate. HepG2 whole cell lysate (ab7900). Paraffin-embedded SAS xenograft
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab97530 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 47 kDa.
IHC-P		1/100 - 1/500.

靶标

功能	Acts as a positive regulator of apoptosis. Also acts as a regulator of cellular reactive oxygen species.
组织特异性	Highly expressed in placenta. Lower levels in heart, lung, skeletal muscle, kidney and pancreas.
序列相似性	Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. DAP kinase subfamily. Contains 1 protein kinase domain.
翻译后修饰	Autophosphorylated.
细胞定位	Nucleus.


图片



Anti-STK17A antibody (ab97530) at 1/1000 dilution + Hep G2 whole cell lysate at 30 μ g

Predicted band size: 47 kDa

7.5% SDS PAGE



Immunohistochemical analysis of paraffin-embedded SAS xenograft, using ab97530 at 1/500 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-STK17A antibody (ab97530)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors