

Anti-STK38L/Ndr2 antibody [3563C2a] ab84873

2 References **1 图像**

概述

| | |
|-------|---|
| 产品名称 | Anti-STK38L/Ndr2抗体[3563C2a] |
| 描述 | 小鼠单克隆抗体[3563C2a] to STK38L/Ndr2 |
| 宿主 | Mouse |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Recombinant fragment |
| 免疫原 | Recombinant fragment corresponding to Human STK38L/Ndr2. Database link: Q9Y2H1 |
| 阳性对照 | Immunised recombinant protein |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | <p>pH: 7.40</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: PBS, 1% BSA, 0.03% Potassium phosphate, 0.812% Sodium chloride, 0.1312% Sodium phosphate, 0.0225% Potassium chloride</p> |
| 纯度 | Protein G purified |
| 纯化说明 | Purified from culture supernatant of hybridoma cultured in a medium containing bovine IgG depleted (approximately 95%) fetal bovine serum. Filtered through a 0.22 µm membrane. |
| 克隆 | 单克隆 |
| 克隆编号 | 3563C2a |
| 同种型 | IgG1 |

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab84873于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | Use at an assay dependent concentration. Predicted molecular weight: 55 kDa. |

靶标

| | |
|-------|--|
| 功能 | Involved in the regulation of structural processes in differentiating and mature neuronal cells. |
| 组织特异性 | Ubiquitously expressed with highest levels observed in the thymus. |
| 序列相似性 | Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. Contains 1 AGC-kinase C-terminal domain. Contains 1 protein kinase domain. |
| 翻译后修饰 | Phosphorylated upon DNA damage, probably by ATM or ATR. |
| 细胞定位 | Cytoplasm. Cytoplasm > cytoskeleton. Associated with the actin cytoskeleton. |

图片



Anti-STK38L/Ndr2 antibody [3563C2a] (ab84873) + immunising recombinant protein

Predicted band size: 55 kDa

Observed band size: 54 kDa

Note: actual band size differs from predicted one, because the WB is depicting the recombinant fragment used as immunogen.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors