

SureLight® APC Anti-T7 tag® antibody ab72565

[1 References](#) [1 图像](#)

概述

产品名称	SureLight® APC Anti-T7 tag®抗体
描述	SureLight® APC兔多克隆抗体to T7 tag®
宿主	Rabbit
偶联物	SureLight® APC. Ex: 652nm, Em: 657nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Species independent
免疫原	Synthetic peptide T7 tag conjugated to a carrier molecule
常规说明	Excitation maximum:650 nm Emission maximum:660 nm
	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.40 Preservative: 0.013% Sodium azide Constituents: 0.1% BSA, 0.87% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
纯度	Size exclusion
纯化说明	Size exclusion chromatography as well as some proprietary methods for purifying the final materials werer used.
克隆	多克隆
同种型	IgG

The Abpromise guarantee **Abpromise™**承诺保证使用ab72565于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

靶标

相关性 The T7 tag is an 11 amino acid peptide encoded in the leader sequence of T7 bacteriophage gene10. This gene encodes a T7 major capsid protein whose function is not clear. The T7 tag serves as a tag in many expression vectors including the pET system that is based on the very efficient T7 RNA polymerase expression system. Monoclonal antibodies specific for T7 tag are an important tool for studying expression of recombinant T7-tagged proteins.

图片



Flow cytometry analysis of a recombinant (N-terminus T7-tag) transmembrane cell surface receptor in CHO cells using ab72565 at a 1/20 dilution (yellow). Black - Mouse IgG SureLight® Allophycocyanin isotype control.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors