

# Anti-Syndecan-1 antibody [B-A38], prediluted ab1410

## 2 References

概述	
产品名称	Anti-Syndecan-1抗体[B-A38], prediluted
描述	小鼠单克隆抗体[B-A38] to Syndecan-1, prediluted
宿主	Mouse
经测试应用	适用于: IHC-Fr, IHC-P
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Syndecan-1. U266 cell line derived from Hodgkins disease
阳性对照	Tonsil
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.30 Constituents: PBS, Tissue culture supernatant
克隆	单克隆
克隆编号	B-A38
骨髓瘤	unknown
同种型	IgG1
轻链类型	unknown

**The Abpromise guarantee** **Abpromise™**承诺保证使用ab1410于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		
IHC-P		

**应用说明**

IHC:

This antibody is pre-diluted for use with an ABC system.

For formalin fixed paraffin embedded sections, perform heat mediated antigen retrieval before commencing with IHC.

**靶标****功能**

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

**序列相似性**

Belongs to the syndecan proteoglycan family.

**细胞定位**

Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors