

Anti-Taurine antibody ab9448

1 References

概述

产品名称	Anti-Taurine抗体
描述	兔多克隆抗体to Taurine
宿主	Rabbit
特异性	The antibody is calibrated against a spectrum of antigens to assure hapten selectivity. Fixed tissue cross-reactivity was tested with known targets at the recommended dilution. No measurable glutaraldehyde-fixed tissue cross-reactivity (<1:1000) was detected against L-alanine, gamma-aminobutyrate, agmatine, guanidine, D/L-arginine, L-citrulline, L-cysteine, D/L-glutamate, D/L-glutamine, glutathione, glycine, L-lysine, L-ornithine, L-serine, L-threonine, L-tryptophan, L-tyrosine.
经测试应用	适用于: IHC-Fr, ICC, IHC-FoFr, Immunomicroscopy
免疫原	Chemical/ Small Molecule corresponding to Taurine (Glutaraldehyde).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Thimerosal (merthiolate) Constituents: 1% Goat Serum, 0.268% PBS
纯度	IgG fraction
克隆	多克隆
同种型	IgG

The Abpromise guarantee **Abpromise™**承诺保证使用ab9448于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.
IHC-FoFr		Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.

靶标

相关性

Taurine is conjugated via its amino terminal group with the bile acids chenodeoxycholic acid and cholic acid to form the bile salts sodium taurochenodeoxycholate and sodium taurocholate. The low pKa (1.5) of taurine's sulfonic acid group ensures that this moiety is negatively charged in the pH ranges normally found in the intestinal tract and thus improves the surfactant properties of the cholic acid conjugate. Taurine has also been implicated in a wide array of other physiological phenomena including inhibitory neurotransmission, long-term potentiation in the striatum/hippocampus, membrane stabilization, feedback inhibition of neutrophil/macrophage respiratory bursts, adipose tissue regulation and calcium homeostasis. The evidence for these claims, when compared against that reported for taurine's role in bile acid synthesis and osmoregulation, is relatively poor. Premature born infants who lack the enzymes needed to convert cystathione to cysteine may become deficient in taurine. Thus, taurine is a dietary essential nutrient in these individuals. Many therapeutic applications of taurine have been investigated. Some conditions that taurine might be useful in treating include: cardiovascular diseases, hypercholesterolemia, epilepsy and other seizure disorders, macular degeneration, Alzheimer's disease, hepatic disorders, alcoholism, and cystic fibrosis. Recent studies show that taurine supplements taken by mice on a high-fat diet reduced their overall weight. Studies have yet to be done on the effect of taurine on obesity in humans.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors