

Anti-TCF19/SC1 antibody ab96828

1 图像

概述

产品名称	Anti-TCF19/SC1抗体
描述	兔多克隆抗体to TCF19/SC1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	<p>Synthetic peptide corresponding to Human TCF19/SC1 aa 296-359. (AAB60363)</p> <p>Sequence:</p> <p>VAACPRKRQWPGFSVMAVTSGSMWPVLAAASRLPGRPTSDAQ</p> <p>GAGLAFSL RSTAKAPSDTPAHE</p> <div>  Run BLAST with  Run BLAST with </div>
阳性对照	293T whole cell lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.00</p> <p>Preservative: 0.01% Thimerosal (merthiolate)</p> <p>Constituents: 89.99% PBS, 10% Glycerol (glycerin, glycerine)</p>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab96828于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 37 kDa.

靶标

功能	Potential trans-activating factor that could play an important role in the transcription of genes required for the later stages of cell cycle progression.
序列相似性	Contains 1 FHA domain. Contains 1 PHD-type zinc finger.
发展阶段	Growth regulated.
翻译后修饰	Phosphorylated upon DNA damage, probably by ATM or ATR.
细胞定位	Nucleus.

图片



Anti-TCF19/SC1 antibody (ab96828) at 1/1500 dilution + 293T whole cell lysate at 30 µg

Predicted band size: 37 kDa

10% SDS-PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors