

### Anti-TdT antibody [41A] ab68886

**1 References**   **3 图像**

#### 概述

产品名称	Anti-TdT抗体[41A]
描述	小鼠单克隆抗体[41A] to TdT
宿主	Mouse
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein. This information is proprietary to Abcam and/or its suppliers.
阳性对照	Thymus.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze.
存储溶液	<p>pH: 7.4</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p>
纯度	IgG fraction
纯化说明	This antibody was purified using a AKTA prime.
克隆	单克隆
克隆编号	41A
同种型	IgG2a

#### 应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab68886于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/100.
WB		1/100. Detects a band of approximately 60 kDa (predicted molecular weight: 58 kDa).

靶标


功能	Template-independent DNA polymerase which catalyzes the random addition of deoxynucleoside 5'-triphosphate to the 3'-end of a DNA initiator. One of the in vivo functions of this enzyme is the addition of nucleotides at the junction (N region) of rearranged Ig heavy chain and T-cell receptor gene segments during the maturation of B- and T-cells.
序列相似性	Belongs to the DNA polymerase type-X family. Contains 1 BRCT domain.
细胞定位	Nucleus.

图片




Anti-TdT antibody [41A] (ab68886) at 1/100 dilution + Thymus extracts

**Predicted band size:** 58 kDa  
**Observed band size:** 60 kDa



Ab68886 (1:100) staining human TdT in human thymus by immunohistochemistry using formalin fixed, paraffin embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TdT antibody [41A] (ab68886)



Ab68886 (1:100) staining TdT in human thymus tissue by immunohistochemistry using formalin fixed, paraffin embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TdT antibody [41A] (ab68886)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors