

# Anti-Testosterone 3 CMO antibody [3T16] ab20250

## 概述

产品名称	Anti-Testosterone 3 CMO抗体[3T16]
描述	小鼠单克隆抗体[3T16] to Testosterone 3 CMO
宿主	Mouse
经测试应用	适用于: RIA
种属反应性	与反应: Species independent
免疫原	Chemical/ Small Molecule corresponding to Testosterone 3 CMO conjugated to bovine serum albumin. Testosterone 3 CMO conjugated to BSA.
常规说明	<p>This antibody can be used for detecting and measuring testosterone levels in serum and tissue.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
Primary antibody说明	This antibody can be used for detecting and measuring testosterone levels in serum and tissue.
克隆	单克隆
克隆编号	3T16
骨髓瘤	P3-NS1/1-Ag4-1
同种型	IgG1

## 应用

## The Abpromise guarantee

**Abpromise™**承诺保证使用ab20250于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
RIA		Use at an assay dependent dilution.

## 靶标

## 相关性

Testosterone-3-carboxy-methyl-oxyme (T-3-CMO). Testosterone belongs to a class of hormones called androgens. These are male hormones responsible for the development of the male reproductive system and secondary male sexual characteristics such as voice depth and facial hair. Testosterone is normally produced by the testes in large quantities in men. It also occurs normally in smaller quantities in women. Several man-made derivatives of testosterone are used to treat advanced disseminated breast cancer in women, especially when cancer has spread to the bones.

## 细胞定位

Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors