

# Anti-TGF beta Receptor III/TGFBR3 antibody ab97459

★★★★★ [1 Abreviews](#) [4 References](#) [1 图像](#)

### 概述

产品名称	Anti-TGF beta Receptor III/TGFBR3抗体
描述	兔多克隆抗体to TGF beta Receptor III/TGFBR3
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Pig 
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
阳性对照	NT2D1, IMR32, U87MG and MCF7 cell lysates
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 78.99% PBS, 1% BSA, 20% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab97459于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	1/500 - 1/3000. Predicted molecular weight: 93 kDa.

## 靶标

功能	Binds to TGF-beta. Could be involved in capturing and retaining TGF-beta for presentation to the signaling receptors.
序列相似性	Contains 1 ZP domain.
翻译后修饰	Extensively modified by glycosaminoglycan groups (GAG).
细胞定位	Secreted. Secreted > extracellular space. Cell membrane. Exists both as a membrane-bound form and as soluble form in serum and in the extracellular matrix.

## 图片



Anti-TGF beta Receptor III/TGFBR3 antibody (ab97459) at 1/1000 dilution + IMR32 whole cell lysate at 30 µg

**Predicted band size:** 93 kDa

7.5% SDS PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors