

# Anti-TRIM30/Rpt1 antibody ab76972

[1 References](#) [2 图像](#)

### 概述

---

|              |   |
|--------------|---|
| <b>产品名称</b>  | Anti-TRIM30/Rpt1 抗体   |
| <b>描述</b>    | 兔多克隆抗体 to TRIM30/Rpt1   |
| <b>宿主</b>    | Rabbit  |
| <b>经测试应用</b> | <b>适用于:</b> WB, IHC-P   |
| <b>种属反应性</b> | <b>与反应:</b> Mouse   |
| <b>免疫原</b>   | Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.   |
| <b>阳性对照</b>  | WB: Mouse spleen tissue lysate. IHC-P: Mouse ovary tissue.  |
| <b>常规说明</b>  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

### 性能

---

|             |   |
|-------------|---|
| <b>形式</b>   | Liquid  |
| <b>存放说明</b> | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| <b>存储溶液</b> | pH: 7.2<br>Preservative: 0.02% Sodium azide<br>Constituent: PBS   |
| <b>纯度</b>   | Immunogen affinity purified   |
| <b>纯化说明</b> | Affinity chromatography purified via peptide column.  |
| <b>克隆</b>   | 多克隆   |
| <b>同种型</b>  | IgG   |

### 应用

---

## The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab76972 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


| 应用    | Ab评论 | 说明  |
|-------|------|---|
| WB    |      | Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 57 kDa. |
| IHC-P |      | Use a concentration of 10 µg/ml.  |

## 靶标

**相关性** Trans acting factor that regulates gene expression of interleukin 2 receptor alpha chain and of human immunodeficiency virus type 1. May affect IL2R alpha and HIV 1 expression through cis acting negative regulatory elements or through competition with proteins that bind to enhancer or activator sequences.

**细胞定位** Nuclear

## 图片




Western blot - Anti-TRIM30/Rpt1 antibody (ab76972)

Anti-TRIM30/Rpt1 antibody (ab76972) at 1 µg/ml + mouse spleen tissue lysate

**Predicted band size:** 57 kDa

**Observed band size:** 57 kDa

**Additional bands at:** 65 kDa, 70 kDa, 80 kDa. We are unsure as to the identity of these extra bands.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TRIM30/Rpt1 antibody (ab76972)

ab76972 at 10 µg/ml staining TRIM30/Rpt1 in mouse ovary tissue section by Immunohistochemistry (Formalin/ PFA fixed paraffin embedded tissue sections).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors