

# Anti-Trypanosoma brucei phosphoglycerate kinase antibody ab89651

★★★★★ [2 Abreviews](#) [1 图像](#)

### 概述

产品名称	Anti-Trypanosoma brucei phosphoglycerate kinase抗体
描述	兔多克隆抗体 to Trypanosoma brucei phosphoglycerate kinase
宿主	Rabbit
经测试应用	适用于: ICC/IF, WB
种属反应性	与反应: Trypanosoma brucei
免疫原	Synthetic peptide derived from the C terminal domain of Trypanosoma brucei phosphoglycerate kinase.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee** **Abpromise™** 承诺保证使用ab89651于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
WB	★★★★★ (1)	1/200 - 1/2000. Predicted molecular weight: 55 kDa.

## 靶标

## 相关性

Belongs to the phosphoglycerate kinase family.

## 图片



**All lanes :** Anti-Trypanosoma brucei phosphoglycerate kinase antibody (ab89651) at 1/200 dilution (in TBS + 5% milk for 12 hours at 4°C)

**Lane 1 :** Whole cell lysate of Trypanosoma brucei - procyclic form

**Lane 2 :** Whole cell lysate of Trypanosoma brucei - bloodstream

Lysates/proteins at 20000000 cells per lane.

## Secondary

**All lanes :** A biotin-conjugated Goat anti-rabbit IgG polyclonal at 1/200 dilution

**Predicted band size:** 55 kDa

**Observed band size:** 65,75 kDa

**Blocking Step:** 5% Milk for 1 hour at 25°C

**Detection Method:** DAB

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors