


Anti-Tryptophan Hydroxylase/TPH antibody ab32821

★★★★★ 1 Abreviews 1 图像

概述

产品名称	Anti-Tryptophan Hydroxylase/TPH抗体
描述	羊多克隆抗体to Tryptophan Hydroxylase/TPH
宿主	Sheep
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: a wide range of other species 
免疫原	Recombinant full length protein corresponding to Rabbit Tryptophan Hydroxylase/TPH.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.50 Constituents: 0.238% HEPES, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 0.01% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab32821于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000. Detects a band of approximately 55 kDa (predicted molecular weight: 51 kDa).

## 靶标

组织特异性	Isoform 2 seems to be less widely expressed than isoform 1.
通路	Aromatic compound metabolism; serotonin biosynthesis; serotonin from L-tryptophan: step 1/2.
序列相似性	Belongs to the bipterin-dependent aromatic amino acid hydroxylase family. Contains 1 ACT domain.

## 图片



Anti-Tryptophan Hydroxylase/TPH antibody (ab32821) at 1/1000 dilution + Human dorsal Raphe nucleus extract at 50 µg

### Secondary

Anti-Sheep HRP labelled.

**Predicted band size:** 51 kDa

**Observed band size:** 55 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors