

Anti-UHRF1BP1L antibody ab117450

[1 References](#) [1 图像](#)

概述

产品名称	Anti-UHRF1BP1L 抗体
描述	兔多克隆抗体 to UHRF1BP1L
宿主	Rabbit
特异性	At least two isoforms of UHRF1BP1L are known to exist; this antibody will recognize both isoforms. ab117450 is predicted to not cross-react with UHRF1BP1.
经测试应用	适用于: WB
种属反应性	与反应: Mouse
免疫原	Synthetic peptide corresponding to Human UHRF1BP1L (internal sequence). Database link: NP_055869
阳性对照	WB: Mouse brain tissue lysate.
常规说明	Avoid repeated freeze thaw cycles and exposure to prolonged high temperatures. The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab117450 was purified by affinity chromatography via a peptide column.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee


Abpromise™承诺保证使用ab117450于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 164 kDa.

靶标

图片



Anti-UHRF1BP1L antibody (ab117450) at 1 µg/ml + mouse brain tissue lysate at 15 µg

Predicted band size: 164 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors