

Anti-UL42 antibody [2H4] ab19311

6 References

概述

产品名称	Anti-UL42抗体[2H4]
描述	小鼠单克隆抗体[2H4] to UL42
宿主	Mouse
特异性	2H4 will recognise UL42 bound to HSV polymerase.
经测试应用	适用于: IP, ELISA, IHC-Fr, WB, ICC/IF
种属反应性	与反应: Herpes simplex virus
免疫原	Baculovirus-expressed HSV DNA polymerase (POL) and POL/UL42 complex.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
克隆	单克隆
克隆编号	2H4
骨髓瘤	Sp2/0
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab19311 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP		
ELISA		
IHC-Fr		
WB		
ICC/IF		

应用说明

ELISA: Use at an assay dependent dilution.
ICC/IF: Use at an assay dependent dilution PMID 17553899.
IHC-Fr: Use at an assay dependent dilution.
IP: Use at an assay dependent dilution.
WB: Use at an assay dependent dilution. Predicted molecular weight: 13.6 kDa.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性 UL42 is the DNA polymerase accessory protein of herpes simplex virus type 1.
细胞定位 Nuclear

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors