

Anti-Ultraspiracle antibody ab106341

2 References

概述

| | |
|-------|---|
| 产品名称 | Anti-Ultraspiracle抗体 |
| 描述 | 兔多克隆抗体to Ultraspiracle |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: ICC, WB |
| 种属反应性 | 与反应: Drosophila melanogaster |
| 免疫原 | Synthetic peptide corresponding to Ultraspiracle (N terminal). Database link: P49700 |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C long term. |
| 存储溶液 | Constituent: Whole serum |
| 纯度 | Whole antiserum |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab106341于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|-----|------|--|
| ICC | | Use at an assay dependent concentration. |
| WB | | Use at an assay dependent concentration. |

靶标

| | |
|-------|--|
| 功能 | Receptor for ecdysone. May be an important modulator of insect metamorphosis. Plays an important part in embryonic and post-embryonic development. Binds to ecdysone response elements (ECRES) such as in the promoter region of s15 chorion gene. |
| 序列相似性 | Belongs to the nuclear hormone receptor family. NR2 subfamily. Contains 1 nuclear receptor DNA-binding domain. |
| 结构域 | Composed of three domains: a modulating N-terminal domain, a DNA-binding domain and a C-terminal ligand-binding domain. |
| 细胞定位 | Nucleus. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors