

Anti-VENTX antibody ab105352

2 References 1 图像

概述	
产品名称	Anti-VENTX抗体
描述	兔多克隆抗体to VENTX
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse
免疫原	Synthetic peptide, 16 amino acid, within the region of the center of Human VENTX (NP_055283).
阳性对照	WB: Mouse brain tissue lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee      Abpromise™承诺保证使用ab105352于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 28 kDa.

## 靶标

功能	May be involved in ventralization.
组织特异性	Expressed in bone marrow of patients recovering from chemotherapy. Also expressed in an erythroleukemia cell line.
序列相似性	Contains 1 homeobox DNA-binding domain.
细胞定位	Nucleus.

## 图片



**Lane 1 :** Anti-VENTX antibody (ab105352) at 1 µg/ml

**Lane 2 :** Anti-VENTX antibody (ab105352) at 2 µg/ml

**All lanes :** Mouse brain tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 28 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors