

Anti-West Nile Virus (E protein) antibody ab25969

概述

产品名称	Anti-West Nile Virus (E蛋白)抗体
描述	兔多克隆抗体to West Nile Virus (E蛋白)
宿主	Rabbit
经测试应用	适用于: ELISA
种属反应性	与反应: West Nile virus
免疫原	Synthetic peptide, corresponding to 16 N terminal amino acids of West Nile Virus (E protein) (Genbank accession no. NP_776014) (Peptide available as ab39903 .)
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab25969于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use a concentration of 1 µg/ml. This will detect 10 ng of free peptide.

靶标

相关性 West Nile Virus (WNV) is a single-stranded RNA (positive sense) virus of the family Flaviviridae, genus Flavivirus, commonly found in Africa, West Asia and the Middle East. It causes disease that is characterized by flu-like symptoms. In a small number of cases this can develop into West Nile Encephalitis (infection of the brain and spinal cord). The viral envelope consists of envelope E and membrane M proteins. The roles of the E protein include virion assembly, recognition of cell receptors, cell endosomal membrane fusion, agglutination of red blood cells, and induction of immune responses.

细胞定位 Cell Membrane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors