

# Anti-White Spot Syndrome Virus VP664 antibody ab26811

## 1 图像

### 概述

产品名称	Anti-White Spot Syndrome Virus VP664抗体
描述	兔多克隆抗体to White Spot Syndrome Virus VP664
宿主	Rabbit
经测试应用	适用于: ELISA
种属反应性	与反应: White spot syndrome virus
免疫原	Synthetic peptide located on the putative domain of White spot syndrome virus VP664.
阳性对照	synthetic WSSV VP664 peptideantigen.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab26811于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ELISA		1/1000.

## 靶标

## 相关性

White spot syndrome virus [WSSV] is one of the most devastating shrimp pathogens, and it has caused serious damage to the shrimp culture industry. This virus can infect several crustacean species, including shrimp, crab, and crayfish. VP664, is one of the viral structural proteins. The existence of the corresponding 18-kb transcript was confirmed by sequencing analysis of reverse transcription-PCR products, which also showed that vp664 was intron-less.

## 图片



ELISA using ab26811 at varying antibody concentrations and antigen concentration at 1000 ng/mL.

ELISA - Anti-White Spot Syndrome Virus VP664  
antibody (ab26811)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors