

Anti-Wnt3a antibody [3A6] ab81614

★★★★★ [2 Abreviews](#) [15 References](#) [1 图像](#)

概述

产品名称	Anti-Wnt3a抗体[3A6]
描述	小鼠单克隆抗体[3A6] to Wnt3a
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length human Wnt3a (aa19-352) purified from E. coli (NP_149122).
阳性对照	WB: A549 cell lysates, Wnt3a transfected 293T cell lysate .
常规说明	<p>This product was changed from ascites to tissue culture supernatant on 24/01/2019. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Ammonium Sulphate Precipitation
克隆	单克隆
克隆编号	3A6
同种型	IgG2a

kappa

The Abpromise guarantee


“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	 (1)	Use at an assay dependent concentration. Detects a band of approximately 39 kDa (predicted molecular weight: 39 kDa).

靶标

功能	Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube.
组织特异性	Moderately expressed in placenta and at low levels in adult lung, spleen, and prostate.
序列相似性	Belongs to the Wnt family.
细胞定位	Secreted > extracellular space > extracellular matrix.

图片



All lanes : Anti-Wnt3a antibody [3A6] (ab81614) at 1/1000 dilution

Lane 1 : 293T cell lysate

Lane 2 : Wnt3a transfected 293T cell lysate

Lysates/proteins at 40 µg per lane.

Secondary

All lanes : goat anti-mouse secondary antibody conjugated to HRP

Predicted band size: 39 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors