


Anti-Wnt3a antibody ab28472

★★★★★ [6 Abreviews](#) [114 References](#) [2 图像](#)

概述

产品名称	Anti-Wnt3a抗体
描述	兔多克隆抗体to Wnt3a
宿主	Rabbit
特异性	Replenishment batches of our polyclonal antibody, ab28472 are tested in WB. Previous batches were additionally validated in ICC/IF. This application is still expected to work and is covered by our Abpromise guarantee. You may also be interested in our alternative recombinant antibody, ab219412 .
经测试应用	适用于: WB, ICC/IF 不适用于: IHC-P
种属反应性	与反应: Mouse, Human 预测可用于: Chicken, Xenopus laevis, Zebrafish, a wide range of other species 
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers. (Peptide available as ab23327)
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituent: PBS

Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our

纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG


The Abpromise guarantee **Abpromise™承诺保证使用ab28472于以下的经测试应用**

应用	Ab评论	说明
WB	★★★★☆ (1)	Use a concentration of 1 µg/ml. Detects a band of approximately 39 kDa (predicted molecular weight: 39 kDa).Can be blocked with <u>Human Wnt3a peptide (ab23327)</u> .
ICC/IF		Use a concentration of 5 µg/ml.


靶标

细胞定位 Secreted > extracellular space > extracellular matrix.

图片



68kDa.



ab28472 stained human embryonic stem cells differentiated into a heterogeneous mix of cells at 5 ug/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors