

# Anti-Yellow Fever Virus antibody [3576] ab22839

★★★★☆ [1 Abreviews](#) [1 References](#)

### 概述

---

<b>产品名称</b>	Anti-Yellow Fever Virus抗体[3576]
<b>描述</b>	小鼠单克隆抗体[3576] to Yellow Fever Virus
<b>宿主</b>	Mouse
<b>特异性</b>	Does not detect other Flaviviruses.
<b>经测试应用</b>	<b>适用于:</b> ELISA
<b>种属反应性</b>	<b>与反应:</b> Yellow fever virus
<b>免疫原</b>	Tissue, cells or virus corresponding to Yellow fever virus Yellow Fever Virus. Whole Yellow fever virus preparation.
<b>常规说明</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	pH: 7.20 Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
<b>纯度</b>	Protein A purified
<b>克隆</b>	单克隆
<b>克隆编号</b>	3576
<b>同种型</b>	IgG2a

### 应用

---

## The Abpromise guarantee

**Abpromise™**承诺保证使用ab22839于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/20 - 1/200.

## 靶标

**相关性** Yellow Fever Virus is a mosquito-borne species of the genus Flavivirus that causes yellow fever in Central and South America and Africa. Yellow fever (also called black vomit or sometimes The American Plague) is an acute viral disease. It is still an important cause of haemorrhage illness in several African and South American countries despite existence of an effective vaccine. In the past it was a source of several devastating epidemics.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors