

Anti-ZCCHC2 antibody ab88756

2 图像

概述	
产品名称	Anti-ZCCHC2抗体
描述	小鼠多克隆抗体to ZCCHC2
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Rat, Human
免疫原	Fragment, corresponding to amino acids 1-362 of Human ZCCHC2
阳性对照	Rat brain tissue lysate and transfected 293T cell lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 100% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee **Abpromise™**承诺保证使用ab88756于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 126 kDa.

靶标

序列相似性

Contains 1 CCHC-type zinc finger.


图片



Anti-ZCCHC2 antibody (ab88756) at 1 µg/ml + Rat brain tissue lysate at 50 µg

Predicted band size: 126 kDa

Observed band size: 126 kDa



All lanes : Anti-ZCCHC2 antibody (ab88756) at 1 µg/ml

Lane 1 : ZCCHC2 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 126 kDa

Observed band size: 40 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors