

Anti-Zebrafish Basement Membrane From Intestinal Epithelial antibody [5F11] ab71312

1 References

概述

产品名称	抗斑马鱼Basement膜From Intestinal Epithelial抗体[5F11]
描述	小鼠单克隆抗体[5F11] to 斑马鱼Basement膜From Intestinal Epithelial
宿主	Mouse
特异性	The transparency of the juvenile zebrafish and its genetic advantages make it an attractive model for the study of intestinal differentiation and renewal. 5F11 is a marker for the basement membrane of the intestinal epithelium.
经测试应用	适用于: ICC/IF, WB, IP
种属反应性	与反应: Zebrafish
免疫原	Tissue, cells or virus corresponding to Zebrafish Zebrafish Basement Membrane From Intestinal Epithelial. Zebrafish Intestine.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	5F11
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab71312于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		
WB		
IP		

应用说明

ICC/IF: Use at an assay dependent dilution.

IP: Use at an assay dependent dilution.

WB: Use at an assay dependent dilution.

This antibody is also suitable for use in histology.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors