

# Anti-ZIP-10 antibody ab83947

[2 References](#) [2 图像](#)

### 概述

---

<b>产品名称</b>	Anti-ZIP-10抗体
<b>描述</b>	兔多克隆抗体to ZIP-10
<b>宿主</b>	Rabbit
<b>特异性</b>	Does not cross-react with ZXinc transporter SLC39A11 (ZIP11)
<b>经测试应用</b>	<b>适用于:</b> WB, ICC
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	Synthetic peptide corresponding to Human ZIP-10 (internal sequence).
<b>常规说明</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

### 应用

---

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab83947于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 94 kDa.
ICC		Use a concentration of 2.5 µg/ml.

## 靶标


功能	May act as a zinc-influx transporter.
序列相似性	Belongs to the ZIP transporter (TC 2.A.5) family.
细胞定位	Membrane.

## 图片



Immunocytochemistry - Anti-ZIP-10 antibody (ab83947)

ab83947 at 2.5µg/ml staining ZIP-10 in Human spleen by Immunohistochemistry.



Western blot - Anti-ZIP-10 antibody (ab83947)

**Lane 1** : Anti-ZIP-10 antibody (ab83947) at 1 µg/ml

**Lane 2** : Anti-ZIP-10 antibody (ab83947) at 2 µg/ml

**All lanes** : Human spleen tissue lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 94 kDa

**Observed band size:** 94 kDa

**Additional bands at:** 50 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors