

Anti-ZIP-12 antibody ab106570

★★★★★ [1 Abreviews](#) [2 References](#) [1 图像](#)

概述

产品名称	Anti-ZIP-12抗体
描述	兔多克隆抗体to ZIP-12
宿主	Rabbit
特异性	At least four isoforms of ZIP-12 are known to exist. ab106570 is predicted not to cross-react with other SLC39 family members.
经测试应用	适用于: ICC/IF, IHC-P, ICC, WB
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide corresponding to Human ZIP-12 (N terminal). Database link: NP_001138667
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab106570 于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.
WB	★★★★★ (1)	Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 77, 73, 46 kDa.

靶标

功能	Acts as a zinc-influx transporter (Potential). May be partly involved in the outbreak of schizophrenia.
组织特异性	Expressed in brain and eye.
序列相似性	Belongs to the ZIP transporter (TC 2.A.5) family.
细胞定位	Membrane.

图片



Lane 1 : Anti-ZIP-12 antibody (ab106570) at 0.5 µg/ml

Lane 2 : Anti-ZIP-12 antibody (ab106570) at 1 µg/ml

All lanes : HepG2 cell lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 77, 73, 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors