

Active human Tollip peptide ab40023

描述

产品名称	Active人Tollip多肽
生物活性	Blocking the activity of TOLLIP antibody.
无动物成分	No
性质	Synthetic
种属	Human

技术指标

Our **Abpromise guarantee** covers the use of **ab40023** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking - Blocking peptide for Anti-Tollip antibody (<u>ab37155</u>)
形式	Liquid

制备和贮存

稳定性和存储	<p>Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.</p> <p>pH: 7.20</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: PBS, 0.1% BSA</p> <p>This product is an active protein and may elicit a biological response in vivo, handle with caution.</p>
--------	---

常规信息

功能	Component of the signaling pathway of IL-1 and Toll-like receptors. Inhibits cell activation by microbial products. Recruits IRAK1 to the IL-1 receptor complex. Inhibits IRAK1 phosphorylation and kinase activity.
序列相似性	<p>Belongs to the tollip family.</p> <p>Contains 1 C2 domain.</p> <p>Contains 1 CUE domain.</p>
翻译后修饰	Phosphorylated by IRAK1 upon stimulation by IL-1 or microbial products.
细胞定位	Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors