

Active Poly (Glu, Tyr) peptide ab204877

2 References

描述

产品名称	Active Poly (Glu, Tyr)多肽
纯度	> 95 % Affinity purified. Purified by affinity chromatography
无动物成分	No
性质	Synthetic
额外的序列信息	Molecular Weight: 20000 - 50000 Da

技术指标

Our **Abpromise guarantee** covers the use of **ab204877** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies
形式	Liquid
补充说明	Ratio for molar concentration 4:1 Glu, Tyr (amount of each amino acid in the solution).

制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle. pH: 7.5 Constituent: 0.394% Tris HCl This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors