

Cow Troponin T1/TNT peptide ab42256

描述

产品名称	牛Troponin T1/TNT多肽
纯度	> 95 % HPLC.
无动物成分	No
性质	Synthetic
种属	Cow
序列	APPPAEVPEVHEEVH
氨基酸	16 to 31

技术指标

Our **Abpromise guarantee** covers the use of **ab42256** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking
形式	Lyophilized
补充说明	This product was previously labelled as Troponin T1

制备和贮存

稳定性和存储	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. Constituent: 0.1% Acetic acid
复溶	Resuspend in 0,1% acetic acid for required concentration

常规信息

功能	Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.
疾病相关	Nemaline myopathy 5
序列相似性	Belongs to the troponin T family.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors