

ELOVL6/LCE peptide ab114079

描述

产品名称	ELOVL6/LCE多肽
无动物成分	No
性质	Synthetic

技术指标

Our **Abpromise guarantee** covers the use of **ab114079** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking - Blocking peptide for Anti-ELOVL6/LCE antibody (<u>ab69857</u>)
形式	Liquid
补充说明	This product was previously labelled as ELOVL6

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 99% PBS, 0.1% BSA
--------	--

常规信息

功能	Condensing enzyme that catalyzes the synthesis of saturated and monounsaturated fatty acids. Highest activity toward C16:0 acyl-CoAs.
组织特异性	Ubiquitous.
序列相似性	Belongs to the ELO family.
翻译后修饰	N-Glycosylated.
细胞定位	Endoplasmic reticulum membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors