

Human beta Defensin 1 peptide ab42255

描述	
产品名称	人beta Defensin 1多肽
纯度	> 95 % HPLC.
无动物成分	No
性质	Synthetic
种属	Human
序列	DHYNCACMGG
氨基酸	1 to 5
技术指标	
Our Abpromise guarantee covers the use of ab42255 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Blocking
形式	Lyophilized
补充说明	Storage: 2-8 °C (lyophilized); - 20 °C (dissolved)
制备和贮存	
稳定性和存储	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. Constituent: 0.1% Acetic acid
复溶	Resuspend in 0,1% acetic acid for required concentration
常规信息	
功能	Has bactericidal activity.
组织特异性	Plasma.
序列相似性	Belongs to the beta-defensin family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors