

Human Perilipin 3/TIP47 peptide ab51558

描述

产品名称	人Perilipin 3/TIP47多肽
无动物成分	No
性质	Synthetic
种属	Human

技术指标

Our **Abpromise guarantee** covers the use of **ab51558** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking - Blocking peptide for Anti-Perilipin 3/TIP47 antibody (ab47639) Western blot
形式	Liquid
补充说明	Previously labelled as Perilipin 3.

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
--------	--

常规信息

功能	Required for the transport of mannose 6-phosphate receptors (MPR) from endosomes to the trans-Golgi network.
序列相似性	Belongs to the perilipin family.
细胞定位	Cytoplasm. Endosome membrane. Lipid droplet. Membrane associated on endosomes. Detected in the envelope and the core of lipid bodies and in lipid sails.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors