

Native Mouse MMP9 protein ab39309

[4 References](#) [1 图像](#)

描述

产品名称	Native小鼠MMP9蛋白
纯度	= 95 % Ion Exchange Chromatography. Purified from cell culture media of BHK cells, free of its endogenous inhibitor, TIMP1, and other matrix metalloproteinases.
表达系统	BHK cells
蛋白长度	Full length protein
无动物成分	No
性质	Native
种属	Mouse

技术指标

Our [Abpromise guarantee](#) covers the use of **ab39309** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE Western blot
形式	Liquid
补充说明	<p>Mouse-MMP9 degrades denatured collagen (gelatin), and a range of extracellular matrix components in-vivo. Unlike MMP2, MMP9 is not constitutively produced by most normal cells, although the enzyme is often over expressed by transformed cells, tumor cells, and cells treated with the phorbol ester PMA. An endogenous inhibitor to MMP9, TIMP1, is often complexed with the enzyme in-vivo, but has been removed from this preparation.</p> <p>ab39309 is a mix of active and latent enzymes, approximately 25% active.</p> <p>For Gelatin zymography the MMP9 does not need to be pre-activated, the SDS in the sample buffer will activate the enzyme sufficiently. For any other assays, or to get maximum activity the MMP9 will need activation, as is is mostly in latent form. APMA at 1 mM works, in a buffer of PBS with 10 mM CaCl₂, 100 μM ZnCl₂, and 0.05% Brij-35, 37°C for a minimum of 1 hour.</p>

制备和贮存

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

pH: 7.40

Constituents: PBS, 50% Glycerol, 1.45% Sodium chloride

常规信息

功能

May play an essential role in local proteolysis of the extracellular matrix and in leukocyte migration. Could play a role in bone osteoclastic resorption. Cleaves KiSS1 at a Gly-Leu bond. Cleaves type IV and type V collagen into large C-terminal three quarter fragments and shorter N-terminal one quarter fragments. Degrades fibronectin but not laminin or Pz-peptide.

组织特异性

Produced by normal alveolar macrophages and granulocytes.

疾病相关

Intervertebral disc disease
Metaphyseal anadysplasia 2

序列相似性

Belongs to the peptidase M10A family.
Contains 3 fibronectin type-II domains.
Contains 4 hemopexin repeats.

结构域

The conserved cysteine present in the cysteine-switch motif binds the catalytic zinc ion, thus inhibiting the enzyme. The dissociation of the cysteine from the zinc ion upon the activation-peptide release activates the enzyme.


翻译后修饰

Processing of the precursor yields different active forms of 64, 67 and 82 kDa. Sequentially processing by MMP3 yields the 82 kDa matrix metalloproteinase-9.
N- and O-glycosylated.

细胞定位

Secreted, extracellular space, extracellular matrix.

图片



Western blot - Recombinant Mouse MMP9 protein (ab39309)

All lanes : Anti-MMP9 antibody ([ab38898](#)) at 2 µg/ml

Lane 1 : Native human MMP9 protein (Proenzyme, monomer) ([ab157344](#))

Lane 2 : Native Mouse MMP9 protein (ab39309)

Lysates/proteins at 0.1 µg per lane.

Performed under reducing conditions.

This blot was produced using a 4-12% Bis-Tris gel under the MOPS buffer system. The gel was run at 200V for 60 minutes before being transferred onto a nitrocellulose membrane at 30V for 70 minutes.

The membrane was then blocked for an hour before being incubated with anti-MMP9 antibody ([ab38898](#); 2 microgram per mL) overnight at 4°C. Antibody binding was detected using infrared labelled goat anti-rabbit (green) antibody (diluted 1:20000) for 1

hour at room temperature before imaging.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors